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1. SAMPLING RECOMMENDATIONS

All the sampling conditions can be consulted via the LNS analysis booklet <https://lns.lu/carnet-danalyse/> and the <https://lns.lu/telechargement/> page, where forms and/or information documents are available.

Patient data and the analysis request must be carefully filled out and verified (matching the prescription), including sample information such as the date (and time of sampling, if relevant), and the identity of person who performed the sampling.

The prescription can serve as the fanalysis request. In this case, the sample information must be entered directly on the prescription.

Below are the recommended sampling conditions according to the type of biological matrix to be sampled:

1.1 WHOLE BLOOD

Matrix	Materials	Additive
Whole blood	Trace Element tube (Royal Blue, dark blue cap) recommended	K ₂ EDTA Sodium Heparin (NaHep)
	EDTA Tube (purple cap)	K ₂ EDTA
	Heparin tube (green cap)	Heparin

1. Samples should preferably be collected in a **trace element K₂EDTA tube /NaHep**.
2. Invert the tube 8 to 10 times to prevent clotting.
3. Avoid opening the tube to limit contamination of the sample.
4. The sample must be transported in the original collection tube. Minimum volume: 0.5 mL.
5. Storage/transport temperature: Refrigerated at +5°C ±3°C is recommended (Ambient or frozen temperatures may be tolerated, depending on the type of analysis to be conducted (e.g., metals).

1.2 URINE

Matrix	Materials	Additive
Urine	clean plastic container	/

1. Collect the samples in a clean **plastic container**.
2. Transfer an aliquot of the well-mixed sample into a transport tube. The minimum volume required for the analyses is :
 - 0.5 mL for metals determination of,
 - 5 mL for **arsenic** speciation,
 - 10 mL for VOC metabolites determination of.
3. Storage/transport temperature: Refrigerated at +5°C ±3°C is recommended (Ambient or frozen temperatures may be tolerated, depending on the type of analysis to be conducted (e.g., metals).

1.3 SERUM

Matrix	Materials	Additive
Serum	Recommended serum trace element tube for inorganic elements	Clot activator
	Gel-free dry tube	Clot activator/no additive

1. Collect the samples in a **royal blue serum tube/gel-free dry tube**. **Glass tubes must not be used to test for aluminium.**
2. After collection, invert the tube 8 to 10 times to mix the clot activator with the blood.
3. Allow the sample to sit at room temperature for 30 minutes, then centrifuge immediately.
4. Transfer serum to a polypropylene transport tube within **2 hours** of collection. **Minimum volume 0.5 mL for inorganic elements and 1 mL for PFAS.**
5. Storage/transport temperature: Refrigerated at +5°C ±3°C is recommended (Ambient or frozen temperatures may be tolerated, depending on the type of analysis to be conducted (e.g., metals).

1.4 ERYTHROCYTES

Matrix	Materials	Additive
Erythrocytes	Trace element tube (Royal Blue) recommended	K ₂ EDTA
	EDTA Tube	K ₂ EDTA
	Heparin tube (green cap)	Heparin

1. Collect the samples in a **K₂EDTA tube**, ideally free of trace elements (EDTA is the preferred anticoagulant to heparin as the latter is only effective for 24 to 36 hours).
2. Invert the tube 8 to 10 times to prevent clotting.
3. The sample must be transported in the original collection tube. **Minimum volume: 4 mL.**
4. Storage/transport temperature: Refrigerated at +5°C ±3°C is recommended (Ambient or frozen temperatures may be tolerated, depending on the type of analysis to be conducted (e.g., metals).

1.5 PLASMA

Matrix	Materials	Additif
Plasma	Trace element tube (Royal Blue) recommended for inorganic elements	K ₂ EDTA
	EDTA Tube	K ₂ EDTA
	Heparin tube (green cap)	Héparin

1. Collect the samples in an anticoagulant treated tubes.
2. Invert the tube 8 to 10 times to prevent clotting.
3. Allow the sample to sit at room temperature for 30 minutes, then centrifuge immediately.
4. Transfer serum to a polypropylene transport tube within 2 hours of collection. Minimum volume 0.5 mL for inorganic elements and 1 mL for PFAS.
5. Storage/transport temperature: Refrigerated at +5°C ±3°C is recommended (Ambient or frozen temperatures may be tolerated, depending on the type of analysis to be conducted (e.g., metals).